

Pcsk1 Antibody, Biotin conjugated



PACO51793

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in the processing of hormone and other protein precursors at sites comprised of pairs of basic amino acid, residues. Substrates include POMC, renin, enkephalin, dynorphin, somatostatin, insulin and AGRP.

Gene ID:

Pcsk1

Uniprot

P28840

Synonyms:

Neuroendocrine convertase 1 (NEC 1) (EC 3.4.21.93) (Prohormone convertase 1) (Proprotein convertase 1) (PC1), Pcsk1, Bdp Nec-1 Nec1

Immunogen:

Recombinant Rat Neuroendocrine convertase 1 protein (111-752AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A