

# PES1 Antibody, FITC conjugated



PACO51800

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Component of the PeBoW complex, which is required for maturation of 28S and 5.8S ribosomal RNAs and formation of the 60S ribosome.

**Gene ID:**

PES1

**Uniprot**

O00541

**Synonyms:**

Pescadillo homolog, PES1

**Immunogen:**

Recombinant Human Pescadillo homolog protein (378-588AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A