

# PLEK Antibody, FITC conjugated



PACO51816

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Major protein kinase C substrate of platelets.

**Gene ID:**

PLEK

**Uniprot**

P08567

**Synonyms:**

Pleckstrin (Platelet 47 kDa protein) (p47), PLEK, P47

**Immunogen:**

Recombinant Human Pleckstrin protein (38-337AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A