

PLEK Antibody, Biotin conjugated



PACO51817

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Major protein kinase C substrate of platelets.

Gene ID:

PLEK

Uniprot

P08567

Synonyms:

Pleckstrin (Platelet 47 kDa protein) (p47), PLEK, P47

Immunogen:

Recombinant Human Pleckstrin protein (38-337AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A