

PLS3 Antibody, HRP conjugated



PACO51819

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Actin-bundling protein found in intestinal microvilli, hair cell stereocilia, and fibroblast filopodia. May play a role in the regulation of bone development.

Gene ID:

PLS3

Uniprot

P13797

Synonyms:

Plastin-3 (T-plastin), PLS3

Immunogen:

Recombinant Human Plastin-3 protein (310-455AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A