

# PLS3 Antibody, FITC conjugated



PACO51820

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Actin-bundling protein found in intestinal microvilli, hair cell stereocilia, and fibroblast filopodia. May play a role in the regulation of bone development.

**Gene ID:**

PLS3

**Uniprot**

P13797

**Synonyms:**

Plastin-3 (T-plastin), PLS3

**Immunogen:**

Recombinant Human Plastin-3 protein (310-455AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A