

PPP1CA Antibody



PACO51834

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:**Gene ID:**

PPP1CA

Uniprot

P62136

Synonyms:

Serine/threonine-protein phosphatase PP1-alpha catalytic subunit (PP-1A) (EC 3.1.3.16),
PPP1CA, PPP1A


Immunogen:

Recombinant Human Serine/threonine-protein phosphatase PP1-α catalytic
subunit protein (192-330AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO51834 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO51834 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO51834 at dilution of 1:100.