

# PPP1CA Antibody, HRP conjugated



PACO51835

---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	PPP1CA
<b>Source:</b>	<b>Uniprot</b>
Rabbit	P62136
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	Serine/threonine-protein phosphatase PP1-alpha catalytic subunit (PP-1A) (EC 3.1.3.16), PPP1CA, PPP1A
<b>Applications:</b>	<b>Immunogen:</b>
ELISA	Recombinant Human Serine/threonine-protein phosphatase PP1-&alpha; catalytic subunit protein (192-330AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A