

PPP1CA Antibody, Biotin conjugated



PACO51837

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

PPP1CA

Uniprot

P62136

Synonyms:

Serine/threonine-protein phosphatase PP1-alpha catalytic subunit (PP-1A) (EC 3.1.3.16), PPP1CA, PPP1A

Immunogen:

Recombinant Human Serine/threonine-protein phosphatase PP1- α ; catalytic subunit protein (192-330AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A