

PRKD3 Antibody, Biotin conjugated



PACO51869

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Converts transient diacylglycerol (DAG) signals into prolonged physiological effects, downstream of PKC. Involved in resistance to oxidative stress.

Gene ID:

PRKD3

Uniprot

O94806

Synonyms:

Serine/threonine-protein kinase D3 (EC 2.7.11.13) (Protein kinase C nu type) (Protein kinase EPK2) (nPKC-nu), PRKD3, EPK2 PRKCN

Immunogen:

Recombinant Human Serine/threonine-protein kinase D3 protein (326-510AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A