

# PRPS2 Antibody, Biotin conjugated



PACO51877

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Catalyzes the synthesis of phosphoribosylpyrophosphate (PRPP) that is essential for nucleotide synthesis.

**Gene ID:**

PRPS2

**Uniprot**

P11908

**Synonyms:**

Ribose-phosphate pyrophosphokinase 2 (EC 2.7.6.1) (PPRibP) (Phosphoribosyl pyrophosphate synthase II) (PRS-II), PRPS2

**Immunogen:**

Recombinant Human Ribose-phosphate pyrophosphokinase 2 protein (11-106AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A