

RAB27A Antibody, HRP conjugated



PACO51911

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a role in cytotoxic granule exocytosis in lymphocytes. Required for both granule maturation and granule docking and priming at the immunologic synapse.

Gene ID:

RAB27A

Uniprot

P51159

Synonyms:

Ras-related protein Rab-27A (Rab-27) (GTP-binding protein Ram), RAB27A, RAB27

Immunogen:

Recombinant Human Ras-related protein Rab-27A protein (2-221AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A