

# RAB27A Antibody, FITC conjugated



PACO51912

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Plays a role in cytotoxic granule exocytosis in lymphocytes. Required for both granule maturation and granule docking and priming at the immunologic synapse.

**Gene ID:**

RAB27A

**Uniprot**

P51159

**Synonyms:**

Ras-related protein Rab-27A (Rab-27) (GTP-binding protein Ram), RAB27A, RAB27

**Immunogen:**

Recombinant Human Ras-related protein Rab-27A protein (2-221AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A