

RPS4Y1 Antibody, FITC conjugated



PACO51932

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytosol, cytosolic small ribosomal subunit, membrane, nucleoplasm, nucleus, polysome, RNA binding, structural constituent of ribosome, multicellular organism development, nuclear-transcribed mRNA catabolic process, nonsense-mediated decay.

Gene ID:

RPS4Y1

Uniprot

P22090

Synonyms:

40S ribosomal protein S4, Y isoform 1 (Small ribosomal subunit protein eS4), RPS4Y1, RPS4Y

Immunogen:

Recombinant Human 40S ribosomal protein S4, Y isoform 1 protein (197-263AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A