

SLIT3 Antibody, Biotin conjugated



PACO51965

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May act as molecular guidance cue in cellular migration, and function may be mediated by interaction with roundabout homolog receptors.

Gene ID:

SLIT3

Uniprot

O75094

Synonyms:

Slit homolog 3 protein (Slit-3) (Multiple epidermal growth factor-like domains protein 5) (Multiple EGF-like domains protein 5), SLIT3, KIAA0814 MEGF5 SLIL2

Immunogen:

Recombinant Human Slit homolog 3 protein (1199-1494AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A