

STAT4 Antibody, FITC conjugated



PACO51984

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Carries out a dual function: signal transduction and activation of transcription. Involved in IL12 signaling.

Gene ID:

STAT4

Uniprot

Q14765

Synonyms:

Signal transducer and activator of transcription 4, STAT4

Immunogen:

Recombinant Human Signal transducer and activator of transcription 4 protein (569-748AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A