

TH Antibody, FITC conjugated



PACO52000

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays an important role in the physiology of adrenergic neurons.

Gene ID:

TH

Uniprot

P07101

Synonyms:

Tyrosine 3-monooxygenase (EC 1.14.16.2) (Tyrosine 3-hydroxylase) (TH), TH, TYH

Immunogen:

Recombinant Human Tyrosine 3-monooxygenase protein (2-197AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A