

TMEM173 Antibody, FITC conjugated



PACO52008

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	TMEM173
Source:	Uniprot
Rabbit	Q86WW6
Isotype:	Synonyms:
IgG	Stimulator of interferon genes protein (hSTING) (Endoplasmic reticulum interferon stimulator) (ERIS) (Mediator of IRF3 activation) (hMITA) (Transmembrane protein 173), TMEM173, ERIS MITA STING
Applications:	Immunogen:
ELISA	Recombinant Human Stimulator of interferon genes protein (207-379AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A