

UGT1A9 Antibody, FITC conjugated



PACO52036

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

UDPGT is of major importance in the conjugation and subsequent elimination of potentially toxic xenobiotics and endogenous compounds. This isoform has specificity for phenols. Isoform 2 lacks transferase activity but acts as a negative regulator of isoform 1.

Gene ID:

UGT1A9

Uniprot

O60656

Synonyms:

UDP-glucuronosyltransferase 1-9 (UDPGT 1-9) (UGT1*9) (UGT1-09) (UGT1.9) (EC 2.4.1.17) (UDP-glucuronosyltransferase 1-I) (UGT-1I) (UGT1I) (UDP-glucuronosyltransferase 1A9) (IugP4), UGT1A9, GNT1 UGT1

Immunogen:

Recombinant Human UDP-glucuronosyltransferase 1-9 protein (277-473AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A