

VEGFB Antibody, Biotin conjugated



PACO52045

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Growth factor for endothelial cells. VEGF-B167 binds heparin and neuropilin-1 whereas the binding to neuropilin-1 of VEGF-B186 is regulated by proteolysis.

Gene ID:

VEGFB

Uniprot

P49765

Synonyms:

Vascular endothelial growth factor B (VEGF-B) (VEGF-related factor) (VRF), VEGFB, VRF

Immunogen:

Recombinant Human Vascular endothelial growth factor B protein (22-207AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A