

WISP3 Antibody, HRP conjugated



PACO52051

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Appears to be required for normal postnatal skeletal growth and cartilage homeostasis.

Gene ID:

WISP3

Uniprot

O95389

Synonyms:

WNT1-inducible-signaling pathway protein 3 (WISP-3) (CCN family member 6), WISP3, CCN6

Immunogen:

Recombinant Human WNT1-inducible-signaling pathway protein 3 protein (201-372AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A