

MLN Antibody, HRP conjugated



PACO52083

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays an important role in the regulation of interdigestive gastrointestinal motility and indirectly causes rhythmic contraction of duodenal and colonic smooth muscle.

Gene ID:

MLN

Uniprot

P12872

Synonyms:

Promotilin [Cleaved into: Motilin; Motilin-associated peptide (MAP)], MLN

Immunogen:

Recombinant Human Promotilin protein (26-115AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A