

# HLA-DQB1 Antibody, Biotin conjugated



PACO52097

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

HLA-DQB1

**Uniprot**

P01920

**Synonyms:**

HLA class II histocompatibility antigen, DQ beta 1 chain (MHC class II antigen DQB1), HLA-DQB1, HLA-DQB

**Immunogen:**

Recombinant Human HLA class II histocompatibility antigen, DQ &beta; 1 chain protein (128-175AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A