

# CHRNB1 Antibody, Biotin conjugated



PACO52101

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

After binding acetylcholine, the AChR responds by an extensive change in conformation that affects all subunits and leads to opening of an ion-conducting channel across the plasma membrane.

**Gene ID:**

CHRNB1

**Uniprot**

P11230

**Synonyms:**

Acetylcholine receptor subunit beta, CHRNB1, ACHRB CHRNB

**Immunogen:**

Recombinant Human Acetylcholine receptor subunit &beta; protein (24-244AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A