

# sox10 Antibody, Biotin conjugated



PACO52113

---

## Product Information

**Size:**

50ug

**Reactivity:**

Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

enteric nervous system development, glial cell differentiation, inner ear development, melanocyte differentiation, negative regulation of neurogenesis, neural crest cell development, neural crest cell migration, olfactory bulb development, oligodendrocyte differentiation, otic vesicle morphogenesis.

**Gene ID:**

sox10

**Uniprot**

Q90XD1

**Synonyms:**

Transcription factor SOX-10, Sox10

**Immunogen:**

Recombinant Zebrafish HMG-box transcription factor Sox10 protein (1-485AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A