

# mitfa Antibody, FITC conjugated



PACO52116

---

## Product Information

**Size:**

50ug

**Reactivity:**

Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

nucleoplasm, RNA polymerase II transcription factor activity, sequence-specific DNA binding, circadian rhythm, melanocyte differentiation, pigmentation, transcription, DNA-templated.

**Gene ID:**

mitfa

**Uniprot**

Q9PWC2

**Synonyms:**

**Immunogen:**

Recombinant Zebrafish Microphthalmia-associated transcription factor a protein (1-412AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A