

mitfa Antibody, Biotin conjugated



PACO52117

Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

nucleoplasm, RNA polymerase II transcription factor activity, sequence-specific DNA binding, circadian rhythm, melanocyte differentiation, pigmentation, transcription, DNA-templated.

Gene ID:

mitfa

Uniprot

Q9PWC2

Synonyms:

Immunogen:

Recombinant Zebrafish Microphthalmia-associated transcription factor a protein (1-412AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A