

**Larval cuticle protein A1A Antibody, FITC conjugated**  
**PACO52132**



---

**Product Information**

**Size:**

50ug

**Reactivity:**

Tenebrio molitor

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Component of the cuticle of the larva of Tenebrio molitor.

**Gene ID:**

**Uniprot**

P80681

**Synonyms:**

Larval cuticle protein A1A (TM-LCP A1A) (TM-A1A)

**Immunogen:**

Recombinant Tenebrio molitor Larval cuticle protein A1A protein (1-174AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A