

PACO52154

---

## Product Information

**Size:**

50ug

**Reactivity:**

Yersinia pestis

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:**

Effector proteins function to alter host cell physiology and promote bacterial survival in host tissues.

**Gene ID:**

yopM

**Uniprot**

P17778

**Synonyms:**

Outer membrane protein YopM, yopM, yop48

**Immunogen:**


Recombinant Yersinia pestis Outer membrane protein YopM protein (1-409AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in Recombinant protein. All lanes: yopM antibody at 2 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. predicted band size: 50 kDa. observed band size: 50 kDa.