

SDC4 Antibody, HRP conjugated



PACO52159

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cell surface proteoglycan that bears heparan sulfate. Regulates exosome biogenesis in concert with SDCBP and PDCD6IP.

Gene ID:

SDC4

Uniprot

P31431

Synonyms:

Syndecan-4 (SYND4) (Amphiglycan) (Ryudocan core protein), SDC4

Immunogen:

Recombinant Human Syndecan-4 protein (19-145AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A