

# SDC4 Antibody, Biotin conjugated



PACO52161

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Cell surface proteoglycan that bears heparan sulfate. Regulates exosome biogenesis in concert with SDCBP and PDCD6IP.

**Gene ID:**

SDC4

**Uniprot**

P31431

**Synonyms:**

Syndecan-4 (SYND4) (Amphiglycan) (Ryudocan core protein), SDC4

**Immunogen:**

Recombinant Human Syndecan-4 protein (19-145AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A