

Cyclomaltodextrin glucanotransferase Antibody



PACO52162

Product Information

Size:

50ug

Reactivity:

Bacillus macerans

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:**Gene ID:****Uniprot**

P31835

Synonyms:

Cyclomaltodextrin glucanotransferase (EC 2.4.1.19) (Cyclodextrin-glycosyltransferase) (CGTase)


Immunogen:

Recombinant Bacillus macerans Cyclomaltodextrin glucanotransferase protein (608-713AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: Cyclomaltodextrin glucanotransferase antibody at 2.5 μ g/ml. Secondary: Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 18 kDa. Observed band size: 18 kDa.