

PACO52170

---

## Product Information

**Size:**

50ug

**Reactivity:**

Escherichia coli

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in the glycolate utilization. Catalyzes the condensation and subsequent hydrolysis of acetyl-coenzyme A (acetyl-CoA) and glyoxylate to form malate and CoA.

**Gene ID:**

glcB

**Uniprot**

P37330

**Synonyms:**

Malate synthase G (MSG) (EC 2.3.3.9), glcB, glc

**Immunogen:**

Recombinant Escherichia coli Malate synthase G protein (2-723AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A