

OSM34 Antibody



PACO52198

Product Information

Size:

50ug

Reactivity:

Arabidopsis thaliana

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

defense response to bacterium, incompatible interaction, defense response to fungus, incompatible interaction, response to salt stress.

Gene ID:

OSM34

Uniprot

P50700

Synonyms:

Osmotin-like protein OSM34, OSM34


Immunogen:

Recombinant Arabidopsis thaliana Osmotin-like protein OSM34 protein (23-244AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: OSM34 antibody at 2 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. predicted band size: 44 kDa. observed band size: 44 kDa.