

dlx2a Antibody



PACO52202

Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the development of the inner ear.

Gene ID:

dlx2a

Uniprot

P50574

Synonyms:

Homeobox protein Dlx2a (DLX-2) (Distal-less homeobox gene 2a), dlx2a, dlx2

Immunogen:

Recombinant Zebrafish Homeobox protein Dlx2a protein (1-270AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A