

petE Antibody, Biotin conjugated



PACO52221

Product Information

Size:

50ug

Reactivity:

Prochlorothrix hollandica

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Participates in electron transfer between P700 and the cytochrome b6-f complex in photosystem I.

Gene ID:

petE

Uniprot

P50057

Synonyms:

Plastocyanin, petE

Immunogen:

Recombinant Prochlorothrix hollandica Plastocyanin protein (35-131AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A