

# Anionic trypsin Antibody, Biotin conjugated



PACO52237

---

## Product Information

**Size:**

50ug

**Reactivity:**

Dog

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

extracellular region, extracellular space, calcium ion binding, serine-type endopeptidase activity, collagen catabolic process, digestion, proteolysis.

**Gene ID:**

**Uniprot**

P06872

**Synonyms:**

Anionic trypsin (EC 3.4.21.4)

**Immunogen:**

Recombinant Dog Anionic trypsin protein (24-247AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A