

htpG Antibody



PACO52238

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Molecular chaperone. Has ATPase activity.

Gene ID:

htpG

Uniprot

P0A6Z3

Synonyms:

Chaperone protein HtpG (Heat shock protein C62.5) (Heat shock protein HtpG) (High temperature protein G), htpG


Immunogen:

Recombinant Escherichia coli Chaperone protein HtpG protein (1-624AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: htpG antibody at 3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. predicted band size: 72 kDa. observed band size: 72 kDa.