

Cationic trypsin Antibody, Biotin conjugated



PACO52253

Product Information

Size:

50ug

Reactivity:

Dog

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

P06871

Synonyms:

Cationic trypsin (EC 3.4.21.4)

Immunogen:

Recombinant Dog Cationic trypsin protein (24-246AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A