

PACO52258

Product Information

Size:

50ug

Reactivity:

Escherichia coli O157: H7

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Channel that opens in response to stretch forces in the membrane lipid bilayer. May participate in the regulation of osmotic pressure changes within the cell.

Gene ID:

mscL

Uniprot

P0A743

Synonyms:

Large-conductance mechanosensitive channel, mscL

Immunogen:

Recombinant Escherichia coli O157: H7 Large-conductance mechanosensitive channel protein (1-136AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A