

# porB Antibody, HRP conjugated



PACO52267

---

## Product Information

**Size:**

50ug

**Reactivity:**

Neisseria meningitidis serogroup B / serotype 15

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Serves as a slightly cation selective porin.

**Gene ID:**

porB

**Uniprot**

E6MZM0

**Synonyms:**

Major outer membrane protein P. IB (PIB) (Protein IB) (Class 3 protein) (Porin), porB

**Immunogen:**

Recombinant Neisseria meningitidis serogroup B / serotype 15 Major outer membrane protein P. IB protein (20-331AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A