

porB Antibody, FITC conjugated



PACO52268

Product Information

Size:

50ug

Reactivity:

Neisseria meningitidis serogroup B / serotype 15

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Serves as a slightly cation selective porin.

Gene ID:

porB

Uniprot

E6MZM0

Synonyms:

Major outer membrane protein P. IB (PIB) (Protein IB) (Class 3 protein) (Porin), porB

Immunogen:

Recombinant Neisseria meningitidis serogroup B / serotype 15 Major outer membrane protein P. IB protein (20-331AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A