

# yvgO Antibody



PACO52270

---

## Product Information

**Size:**

50ug

**Reactivity:**

Bacillus subtilis

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

yvgO

**Uniprot**

O32211

**Synonyms:**

Stress response protein YvgO, yvgO

**Immunogen:**

Recombinant Bacillus subtilis Stress response protein YvgO protein (27-161AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A