

Fgf15 Antibody, FITC conjugated



PACO52276

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in the suppression of bile acid, biosynthesis through down-regulation of CYP7A1 expression.

Gene ID:

Fgf15

Uniprot

O35622

Synonyms:

Fibroblast growth factor 15 (FGF-15), Fgf15

Immunogen:

Recombinant Mouse Fibroblast growth factor 15 protein (26-218AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A