

# Fgf15 Antibody, Biotin conjugated



PACO52277

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Involved in the suppression of bile acid, biosynthesis through down-regulation of CYP7A1 expression.

**Gene ID:**

Fgf15

**Uniprot**

O35622

**Synonyms:**

Fibroblast growth factor 15 (FGF-15), Fgf15

**Immunogen:**

Recombinant Mouse Fibroblast growth factor 15 protein (26-218AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A