

ZNF816 Antibody



PACO52290

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF816

Uniprot

Q0VGE8

Synonyms:

Zinc finger protein 816, ZNF816, ZNF816A


Immunogen:

Recombinant Human Zinc finger protein 816 protein (1-200AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO52290 at dilution of 1:100.



Immunofluorescent analysis of PC-3 cells using PACO52290 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human colon cancer using PACO52290 at dilution of 1:100.