

# OVGP1 Antibody, FITC conjugated



PACO52300

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Binds to oocyte zona pellucida in vivo. May play a role in the fertilization process and/or early embryonic development.

**Gene ID:**

OVGP1

**Uniprot**

Q12889

**Synonyms:**

Oviduct-specific glycoprotein (Estrogen-dependent oviduct protein) (Mucin-9) (Oviductal glycoprotein) (Oviductin), OVGP1, MUC9 OGP

**Immunogen:**

Recombinant Human Oviduct-specific glycoprotein protein (484-678AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A