

OVGP1 Antibody, Biotin conjugated



PACO52301

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds to oocyte zona pellucida in vivo. May play a role in the fertilization process and/or early embryonic development.

Gene ID:

OVGP1

Uniprot

Q12889

Synonyms:

Oviduct-specific glycoprotein (Estrogen-dependent oviduct protein) (Mucin-9) (Oviductal glycoprotein) (Oviductin), OVGP1, MUC9 OGP

Immunogen:

Recombinant Human Oviduct-specific glycoprotein protein (484-678AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A