

NMI Antibody, Biotin conjugated



PACO52325

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in augmenting coactivator protein recruitment to a group of sequence-specific transcription factors. Augments cytokine-mediated STAT transcription. Enhances CBP/p300 coactivator protein recruitment to STAT1 and STAT5.

Gene ID:

NMI

Uniprot

Q13287

Synonyms:

N-myc-interactor (Nmi) (N-myc and STAT interactor), NMI

Immunogen:

Recombinant Human N-myc-interactor protein (62-202AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A