

MYLK Antibody, FITC conjugated



PACO52352

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

MYLK

Uniprot

Q15746

Synonyms:

Myosin light chain kinase, smooth muscle (MLCK) (smMLCK) (EC 2.7.11.18) (Kinase-related protein) (KRP) (Telokin) [Cleaved into: Myosin light chain kinase, smooth muscle, deglutamylated form], MYLK, MLCK MLCK1 MYLK1

Immunogen:

Recombinant Human Myosin light chain kinase, smooth muscle protein (290-398AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A